THUASNE

Ortel

fr	Bande herniaire	5
en	Hernia bandage	6
de	Bruchband	7
nl	Breukband	8
it	Cinto erniario	9
es	Faja herniaria	10
pt	Banda para hérnia	11
da	Brokbind	12
fi	Tukivyö tyrälle	13
sv	Bråckband	14
el	Κηλεπίδεσμος	15
cs	Kýlní pás	16
pl	Orteza przepuklinowa	17
lv	Trūces bandāža	19
lt	Išvaržos tvarstis	20
et	Songavöö	21
sl	Oporni pas za hernijo	22
sk	Prietržový pás	23
hu	Sérvkötő	24
bg	Лента за херния	25
ro	Bandă pentru hernie	26
ru	Грыжевой бандаж	27
hr	Povez za herniju	28
zh	疝带	29
ar	حزام الفتق	31

	75 - 85 cm	1
	86 - 97 cm	2
	98 - 110 cm	3
	111 - 125 cm	4

fr	Tour de taille
en	Waist measurement
de	Taillenumfang
nl	Tailleomvang
it	Circonferenza vita
es	Contorno de cintura
pt	Perímetro de cintura
da	Taljemål
fi	Vyötärön ympärys
sv	Midjemått
el	Περίμετρος μέσης
cs	Obvod pasu
pl	Obwód w talii
lv	Vidukļa apkārtmērs
lt	Juomens apimtis
et	Vööümbermõõt
sl	Obseg pasu
sk	Obvod pása
hu	Derékbőség
bg	Обиколка на талията
ro	Circumferința taliei
ru	Обхват талии
hr	Opseg struka
zh	腰围
ar	محيط الخصر

Soutien
Support
Stütze
Ondersteuning
Sostegno
Sujeción
Suporte
Støtte
Tuki
Stöd
Υποστήριξη
Podpora
Podtrzymywanie
Atbalsts
Palaikymas
Toetus
Podpora
Podpora
Megtámasztás
Поддръжка
Sustinere
Поддержка
Potpora
支撐
الدعم

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HERNIA BANDAGE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Strong elastic fabric. ❶

Removable repositionable pads. ❷

Interchangeable adjustable straps. ❸

Exclusive buckle system to fit right or left side. ❹

Textile components: cotton - elastodiene - polyamide - polyester.

Non-textile components: iron - polyamide - polyurethane - silicone.

Properties/Mode of action

The device ensures an adaptable support to the inguinal area thanks to the removable repositionable pads.

Indications

Unilateral ☺ or bilateral ☹ inguinal hernia.

Pre-operative wear.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Do not use in the event of non-reducible hernia.

Precautions

Reduce the hernia before fitting the truss. To do this, it is recommended that you lie on your back and massage the hernia.

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

Do not wear the product in a medical imaging machine.

Do not wear the product while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended that a healthcare professional supervises the first application.

It is recommended to adequately tighten the device to achieve support without restricting blood circulation.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Atrophy of the spermatic cord.

Fusion of the hernia sac.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Adjustment of the device

Prior instructions for the use of the unilateral truss :

The unilateral truss may be used for either a right-side hernia or a left-side hernia.

• Right-side hernia :

The unilateral truss is supplied in the "right-side hernia" position.

Simply pull back the two plastic guides, one after the other.

• Left-side hernia :

Remove strap, then position it on the other side.

Using the two plastic guides, pull off buckle, turn it round and replace it (metal clip towards the outside).

Pull back the two plastic guides one after the other. Turn round the fastener.

The truss is in the "left-side hernia" position.

Fitting the device

Place the pad(s) on the inside surface, against the hernia(s). ❹

Adjust the truss to your waist size by moving fastener. ☹ The truss should be firmly fitted, without being too tight.

Adjust the length of the thigh straps. ☹

Check that the pads are properly positioned.

Care

Hand wash. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products

containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Egyoldali **©** vagy kétoldali **©** lágyéksérv.