THUASNE

Cemen

fr	Bande-ceinture thoracique et abdominale	6
en	Thoracic and abdominal belt-bandage	7
de	Brust- und Bauchgurtband	9
nl	Gordelband voor borstkas en buik	10
it	Fascia-cintura toracica e addominale	12
es	Faja-cinturón torácica y abdominal	13
pt	Banda-cinta torácica e abdominal	15
da	Bryst- og mavebælte	16
fi	Tukiside ylävartalolle ja vatsalle	18
sv	Bältesband för bröstkorgen och buken	19
el	Ταινία-ζώνη θωρακικής και κοιλιακής χώρας	21
cs	Hrudní a břišní pás	22
pl	Orteza piersiowa i brzuszna	24
lv	Krūškurvja un vēdera bandāža-josta	25
lt	Torakalinis ir abdominalinis tvarstis-diržas	27
et	Rindkere ja kõhu vöö	28
sl	Oporni pas za trebušni predel in prsni koš	30
sk	Hrudný a brušný pás	31
hu	Mellkasi és hasi öv	33
bg	Лентовиден колан за гърдите и коремната област	34
ro	Centură toracică și abdominală de tip bandă	36
ru	Бандаж-корсет для груднопоясничного и брюшного отдела	37
hr	Torakalno-abdominalni pojas/povez	39
zh	胸腹支撐帶	40
ar	حزام-ضمادة للصدر والبطن	43


fr	Compression	Soutien
en	Compression	Support
de	Kompression	Stütze
nl	Compressie	Ondersteuning
it	Compressione	Sostegno
es	Compresión	Sujeción
pt	Compressão	Suporte
da	Kompression	Støtte
fi	Kompressio	Tuki
sv	Kompression	Stöd
el	Συμπίεση	Υποστήριξη
cs	Komprese	Podpora
pl	Kompresja	Podtrzymywanie
lv	Kompresija	Atbalsts
lt	Kompresija	Palaikymas
et	Kompressioon	Toetus
sl	Kompresija	Podpora
sk	Kompresia	Podpora
hu	Szorítás	Megtámasztás
bg	Компресия	Поддръжка
ro	Compresie	Susținere
ru	Компрессия	Поддержка
hr	Kompresija	Potpورا
zh	压紧	支撑
ar	الضغط	الدعم

18 cm

25 cm

1,25 m - 2,5 m

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Longueurs : 1,25 m - 2,5 m

Heights: 18 cm - 25 cm

Lengths: 1.25 m - 2.5 m

Höhen: 18 cm - 25 cm

Längen: 1,25 m - 2,5 m

Hoogtes: 18 cm - 25 cm

Lengtes: 1,25 m - 2,5 m

Altezze: 18 cm - 25 cm

Lunghezze: 1,25 m - 2,5 m

Alturas: 18 cm - 25 cm

Longitudes: 1,25 m - 2,5 m

Alturas: 18 cm - 25 cm

Comprimentos: 1,25 m - 2,5 m

Højde: 18 cm - 25 cm

Længder: 1,25 m - 2,5 m

Korkeudet: 18 cm - 25 cm

Pituudet: 1,25 m - 2,5 m

Höjd: 18 cm - 25 cm

Längder: 1,25 m - 2,5 m

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Výšky: 18 cm - 25 cm

Délka: 1,25 m - 2,5 m

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Długości: 1,25 m - 2,5 m

Augstumi: 18 cm - 25 cm

Garumi: 1,25 m - 2,5 m

Aukštis: 18 cm - 25 cm

Ilgiai: 1,25 m - 2,5 m

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Pikkused: 1,25 m - 2,5 m

Višina: 18 cm - 25 cm

Dolžina: 1,25 m - 2,5 m

Výšky: 18 cm - 25 cm

Dĺžky: 1,25 m - 2,5 m

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Hossz: 1,25 m - 2,5 m

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Conserver cette notice.

en

THORACIC AND ABDOMINAL BELT-BANDAGE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Elastic fabric. ①

Thumb-loops. ②

Hooks fastening. ③

Textile components: elastane - polyamide - polyester - viscose.

Non-textile components: stainless steel.

Properties/Mode of action

The device ensures strong support and adaptable compression to the abdomen and/or the rib cage thanks to the elastic fabric.

Indications

Rib fracture.

Postoperative complications (i.e. hysterectomy, bariatric surgery, spinal surgery, cesarean, laparotomy).

Tetraplegic spinal cord injury.

Orthostatic hypotension.

Respiratory rehabilitation.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Precautions

Warning, The very sharp hooks must be handled with care to prevent any risk of injury.

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

Do not wear the product in a medical imaging machine.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional. For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support without excessive compression.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

In postoperative conditions, the restrictive effect of the device could alter respiratory function and decrease lung capacity and tidal volume.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Unroll a section of bandage approximately 80 cm long.

Place the end of the bandage on the front of the body, wind round once, stretching the bandage to obtain the required compression. Ⓐ The label must be positioned outside.

Continue winding, always keeping the bandage stretched to the same level. Ⓑ

Fasten the end of the bandage by attaching the hooks directly onto the fabric using the thumb-loops. Ⓒ

After use, re-roll the bandage, with the side bearing the hooks towards the outside.

Care

Machine washable at 40 °C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Ironing with a cool iron possible. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. Once dry, roll the bandage with the hooks towards the outside.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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وإصل بسط الحزام مع الحفاظ على نفس قوة الجذب. ©
Thuasne SAS – SIREN/RCS Nanterre 542 091 186